



Getting your OxyGo FIT™ ... IS EASY!

Now with the NEW "My OxyGo" APP!



Accessory Options

available through your provider



- 1 Rechargeable Lithium-Ion Battery
Single: #1400-2010-4 Double: #1400-2010-8
- 2 Desktop Battery Charger: #1400-2030
- 3 Car/Boat/RV DC Power Cable: #1400-1050
- 4 Carry Strap (only): #1170-2415
- 5 Wall Charger AC Power Supply: #1400-2040
- 6 Protective Bag: #1170-2410

Your OxyGo FIT Provider:

Empty box for provider information.

MADE IN THE USA

OxyGo is made in the U.S.A. Printed in the U.S.A.

PN B102-0219



Keep Going™

Introducing the newest member of the OxyGo® family...

OxyGo FIT™

The 3 setting Intelligent Pulse dose OxyGo FIT is the smallest and lightest member of the OxyGo family! OxyGo FIT is about half the size of OxyGo, and is small enough and quiet enough to be worn on a patient's hip. Weighing in at 2.8 lbs.* with up to a 5 hour battery life** the OxyGo FIT Keeps Going—Everywhere.

* Single battery ** Double battery



Take it with you—
anywhere!

- Extraordinarily Compact and Quiet
- Charge and Use Anywhere... 24/7!
- 3 Pulse Dose Flow Settings
- Bluetooth-enabled units now connect to your mobile device with the **My OxyGo App***



OxyGo FIT Comes With:

- A durable, black protective bag
- Carry strap
- AC charger for home, DC charger for your car
- Easy to follow instructions!

PLUS – The NEW “My OxyGo” Mobile App!

Check your Bluetooth-enabled OxyGo FIT's battery life, current settings, upcoming maintenance, FAQs — and MORE with our new mobile app!

All the information you need to stay mobile and active is now in the palm of your hand. With the My OxyGo App, you also receive important updates and travel tips for your OxyGo FIT — 24/7!



DOWNLOAD IT TODAY!



Size, Weight	With Standard Battery: (L) 5.91" x (W) 2.68" x (H) 7.2" • 2.8 lbs. / 1.27 kg. With Double Battery: (L) 5.91" x (W) 2.68" x (H) 7.79" • 3.29 lbs. / 1.49 kg.
Oxygen Flow	Per setting: 1 (210ml/min), 2 (420ml/min), 3 (630ml/min)
Power	AC Power Supply: 100-240 V, 50-60 Hz (auto-sensing to allow worldwide use), DC Power for mobile use in car
Battery	Single Battery: Up to 2.7 hours Recharge up to 3 hours with AC or DC power Double Battery: Up to 5 hours Recharge up to 5 hours with AC or DC power
Noise	40 Decibels (at flow setting 2)
Warranty	3 year lifetime limited**
Operation	Simple control functions and easy-to-read LCD display
Use	Designed for 24/7 use, at home or away!

* Bluetooth capability available on OxyGo FIT POCs produced after 2/19

** See OxyGo Warranty Statement for details